
hildren and adolescent sizes

2013-2014

Important
Before printing, turn off
"fit to page" in your print
options dialogue box.*



Length Measurement

- 1
- Place bottom edge of this sheet against wall on a hard surface.
- Tape down to prevent any mis-sizing



- Place heel of AFO on sizing chart **against wall**
- Be certain to include posterior stop
- If foot is longer than brace, include foot




- Trace AFOs for accuracy
- Match longest point to number on right
- If you are between or slightly above line, round up

4 IMPORTANT


Necessary Adjustments for Circumference For all AFO's, TCFlex, DAFO

For all AFO's, TCFlex, DAFO's, KAFO's, SMO's with:

- Solid ankle, instep strap, articulating hinges, posterior stop *Increase 1 full size*
- Turbo boot, postrior stop, hinges, footplates, instep, forefoot and toe strapping: *Increase 2 sizes*



PLACE HEEL OF AFO AGAINST WALL



KEEPINGPAC Children's Orthopedic Footw

Keeping Pace, Inc.

contact us: footwear@keepingpace.com P (888) 526-0020 F (978) 526 0025 USA

> visit us on the web: www.keepingpace.com